



ecore™

Sustainable Solutions for Health, Productivity and the Environment.



Enviro Care® Neutral Disinfectant



Enviro Care[®]

Neutral Disinfectant

Where to Use

Recommended for use in hospitals, nursing homes, schools, colleges, commercial and industrial institutions, veterinary clinics, sports, fitness, and recreation applications. Disinfects, cleans and deodorizes floors, walls, metal surfaces, stainless steel surfaces, glazed porcelain, plastic surfaces and other hard non-porous surfaces.

Benefits

- Hard Water Tolerance
- NPE Free
- Organic Soil Tolerance
- Neutral formula will not dull floor finishes
- Cleans, Disinfects, Deodorizes
- Pseudomonicidal, Fungicidal, Virucidal

Dilution Ratios

Dilute: (1:64) 2 oz per gallon/60 ml per 3.8 litres of water.

Directions for Use

Application: Remove heavy soil deposits from surface. Then thoroughly wet surface. The solution can be applied with a cloth, mop, sponge, or coarse spray, or soaking. Let solution remain on surface for a minimum of 10 minutes. Rinse or allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.

For Human Coronavirus: treated surfaces must remain wet for 1 minute. Wipe dry with a clean cloth, sponge or mop or allow to air dry.

For Influenza Virus Type A: treated surfaces must remain wet for 2 minutes. Wipe dry with a clean cloth, sponge or mop or allow to air dry.

Kills HIV-1 (AIDS virus) on pre-cleaned, environmental surfaces or objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human immunodeficiency virus Type I (HIV-1) (Associated with AIDS).

Effective Against:

- HIV-1 (AIDS virus).
- Hepatitis B (HBV).
- Hepatitis C Virus (HCV).
- Staphylococcus aureus - Methicillin resistant strain (MRSA).
- Escherichia coli (antibiotic resistant).
- Aspergillus Niger, Candida albicans, Trichophyton mentagrophytes (a cause of ringworm).
- Pseudomonas aeruginosa (antibiotic resistant)
- Streptococcus faecalis (antibiotic resistant)
- Klebsiella pneumoniae (antibiotic resistant)

Properties	
Color	Blue
Fragrance	Mild
Viscosity	Water, thin
VOC's (% by weight) EPA Method 24	0.5%
pH (Concentrated)	8.2
Readily Biodegradable	Yes
Phosphate Free	Yes

Available in 1 gallon
containers from Ecore.

1 gallon of Enviro Care makes 64 gallons
of disinfectant.

For handling and precautionary information, reference
the product label and safety data sheet.

For more information about
RMC call 1.800.836.1633 or go
to www.rochestermidland.com

Enviro Care[®]

Neutral Disinfectant

Efficacy Data

Bactericidal

Enviro Care Neutral Disinfectant is bactericidal against many gram-negative pathogens according to the use of the dilution tests AOAC, modified in the presence of serum and hardness of organic water 5% 400 ppm (calculated in CaCO₃).

Dilution 1:64 (2 oz per gallon / 60 mL for 3.8 liters of water).

Acinetobacter baumannii
Acinetobacter calcoaceticus
Bordetella bronchiseptica
Chlamydia psittaci
Enterobacter aerogenes
Enterobacter cloacae
Enterobacter cloacae -NDM-1
Enterococcus faecalis (vancomycin resistant-VRE)
Escherichia coli
Escherichia coli-NDM-1
Fusobacterium necrophorum
Klebsiella pneumoniae
Klebsiella pneumoniae-NDM-1
Legionella pneumophila
Listeria monocytogenes
Pasteurella multocida
Proteus mirabilis
Proteus vulgaris
Pseudomonas aeruginosa
Salmonella choleraesuis
Salmonella enterica
Salmonella enteritidis
Salmonella typhi
Salmonella typhimurium
Serratia marcescens
Shigella flexneri
Shigella sonnei

Community-associated methicillin-resistant Staphylococcus aureus - (CA-MRSA) (NRS384) (USA300)
Community-associated methicillin-resistant Staphylococcus aureus - (CA-MRSA) (NRS123) (USA400)
Staphylococcus aureus-résistant à la méthicilline (MRSA)
Staphylococcus epidermidis
Streptococcus pyogenes (Strep)
Staphylococcus aureus-Vancomycin Intermediate Resistant (VISA)
Streptococcus faecalis

Virucidal

Enviro Care Neutral Disinfectant is virucidal against the following viruses after having been tested according to the virucide classification, according to the use of the dilution tests CaCO₃ modified in the presence of serum and hardness of the organic water 5% 400 ppm (calculated in CaCO₃).

Dilution 1:64 (2 oz per gallon / 60 mL in 3.8 liters of water).

Adenovirus Type 7 - 8 oz. Per gallon
Adenovirus Type 4
Avian Influenza (H5N1)
Avian polyomavirus
Canine distemper

Coronavirus (SARS)- (cause of acute respiratory syndrome)
Equine Herpes Virus Type 1
Feline coronavirus
Feline leukemia
Feline panleukopenia virus
Feline picornavirus
Gastroenteritis transmissible (TGE)
Hepatitis B Virus
Hepatitis C virus
Herpes Simplex Type I
Herpes Simplex Type II
HIV-1 (AIDS)
Human coronavirus
Infectious bovine rhinotracheitis
Infectious bronchitis (Avian IBV)
Influenza A / Hong Kong (H1N1)-Swine flu
Newcastle disease
Porcine parvovirus
Pseudorabies (PRV)
Rabies
Respiratory syncytial virus (RSV)
Rotavirus [WA]
Rubella virus
SARS Associated Coronavirus
Vaccinia (Pox virus)

Fungicidal

Enviro Care Neutral Disinfectant is effective against the following fungi according to the AOAC fungicide test modified in the presence of 5% iron (BIOSTRESS load) and water hardness of 400 ppm (CaCO₃ calculation) at 20°C.

Dilution 1:64 (2 oz per gallon / 60 ml in 3.8 liters of water).

Aspergillus niger
Trichophyton mentagrophytes (athlete's foot)
Candida albicans

Sanitizing

Enviro Care Neutral Disinfectant is a sanitizer for hard nonporous non-food contact surfaces according to the method for Sanitizers - Non-Food Contact Surfaces, modified in the presence of 400 ppm hard water plus 5% organic serum in 60 seconds against:

Enterobacter aerogenes
Staphylococcus aureus

TO SANITIZE NON-FOOD CONTACT HARD NON-POROUS SURFACES: Add 2 oz. of this product per gallon of water to sanitize hard nonporous surfaces. Treated surfaces must remain wet for 60 seconds. Wipe dry with a sponge, mop, or cloth or allow to air dry. Do not use in food preparation, food handling and food storage areas. This product is not to be used on glassware, dishware, cookware or utensils.